

## Fenoprofen Calcium

Catalog No: #S6497


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	Fenoprofen Calcium
Brief Description	Inhibitors
Purification	/
Target Name	Others
Calculated MW	522.6
Formulation	C30H26CaO6
Storage	3 years -20°C powder; 2 years -80°C in solvent;

## Images



## Product Description

Research Area: Other SMILES: [Ca++].CC(C([O-])=O)C1=CC(OC2=CC=CC=C2)=CC=C1.CC(C([O-])=O)C1=CC(OC2=CC=CC=C2)=CC=C1  
 |c:11,13,15,17,30,32,34,36,t:5,9,24,28|Pathways: Others Receptor: Others Boiling pt: Melting pt: Solubility: DMSO: 105 mg/mL (200.91 mM) Ethanol: 10 mg/mL (19.13 mM) Water:

Note: This product is for in vitro research use only