

Fenticonazole Nitrate

Catalog No: #S6499


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Fenticonazole Nitrate
Brief Description	Inhibitors
Purification	98.00%
Target Name	Others
Calculated MW	518.41
Formulation	C ₂₄ H ₂₀ Cl ₂ N ₂ O ₃ HNO ₃
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research

Area: Immunology SMILES: ClC1=CC(Cl)=C(C=C1)C(CN1C=CN=C1)OCC1=CC=C(SC2=CC=CC=C2)C=C1 Pathways: Others Receptor: Others Boiling pt: Melting pt: Solubility: Ethanol: 2 mg/mL (3.85 mM); DMSO: 104 mg/mL (200.61 mM) Appearance: Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several months.

References

1. Fernandez-Alba J, et al. J Chemother. 2004 Apr; 16(2):179-86.

Note: This product is for in vitro research use only