

Fluocinonide

Catalog No: #S6502


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Fluocinonide
Brief Description	Inhibitors
Purification	99.00%
Target Name	GR inhibitor
Calculated MW	494.52
Formulation	C26H32F2O7
Storage	3 years -20°C powder;2 years -80°C in solvent;

Images



Product Description

Research Area:Endocrine

systemSMILES:CC(=O)OCC(=O)[C@@]12[C@@H](C[C@@H]3[C@@]1(C[C@@H]([C@]4([C@H]3C[C@@H](C5=CC(=O)C=C[C@@]54C)F)F)O)C)OC(O2)(C)CPathways:MetabolismReceptor:Glucocorticoid ReceptorBoiling pt:Melting pt:Solubility:DMSO: 99 mg/mL (200.19 mM)Appearance:Remark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

References

1. Gans EH, et al. J Drugs Dermatol, 2008, 7(1), 28-32.

Note: This product is for in vitro research use only