

MRP-S34 Polyclonal Antibody

Catalog No: #41182

Package Size: #41182-1 50ul #41182-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	MRP-S34 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB ELISA
Species Reactivity	Hu Ms
Specificity	MRP-S34 Polyclonal Antibody detects endogenous levels of MRP-S34 protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the C-terminal region of human MRP-S34.
Target Name	MRP-S34
Other Names	MRPS34; 28S ribosomal protein S34; mitochondrial; MRP-S34; S34mt
Accession No.	Swiss-Prot: P82930NCBI Gene ID: 65993
Uniprot	P82930
GeneID	65993;
SDS-PAGE MW	26kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year


Application Details

Western Blot: 1/500 - 1/2000.


ELISA: 1/10000.

Not yet tested in other applications.

Images



Western Blot analysis of Jurkat cells using MRP-S34 Polyclonal Antibody



Western Blot analysis of Jurkat HuvEc cells using MRP-S34 Polyclonal Antibody

Note: This product is for in vitro research use only