

Halobetasol Propionate

Catalog No: #S6529


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Halobetasol Propionate
Brief Description	Inhibitors
Purification	99.00%
Target Name	Phospholipase inhibitor
Calculated MW	484.96
Formulation	C25H31ClF2O5
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Metabolism

systemSMILES:[H][C@@]12C[C@H](C)[C@](OC(=O)CC)(C(=O)CC)[C@@]1(C)C[C@H](O)[C@@]1(F)[C@@]2([H])C[C@H](F)C2=CC(=O)C=C[C@@]12C |r,c:33,t:29|Pathways:MetabolismReceptor:PLA2Boiling pt:Melting pt:Solubility:DMSO:97 mg/mL (200.01 mM)Ethanol:38 mg/mL (78.35 mM)Water:

References

1. Mohandas S, et al. Indian J Dermatol Venereol Leprol, 2009, 75(2), 186-187.

Note: This product is for in vitro research use only