

Heparin sodium

Catalog No: #S6531


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Heparin sodium
Brief Description	Inhibitors
Purification	98.00%
Target Name	Thrombin inhibitor
Calculated MW	6000-20000
Formulation	(C ₁₂ H ₁₆ NS ₂ Na ₃) ₂₀
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research

Area: Immunology SMILES: C([C@@H]1C([C@@H](C(C(O1)O)NS(=O)(=O)[O-])O)O[C@H]2C([C@H](C(C(O2)C(=O)[O-])O)O)OS(=O)(=O)[O-])OS(=O)(=O)[O-])Pathways: Others Receptor: Thrombin III Boling pt: Melting pt: Solubility: DMSO:

References

1. Capila I, et al. Angew Chem Int Ed Engl. 2002, 41(3), 391-412.

Note: This product is for in vitro research use only