

## Hesperadin

Catalog No: #S6532


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	Hesperadin
Brief Description	Inhibitors
Purification	98.00%
Target Name	Aurora Kinase inhibitor
Calculated MW	516.65
Formulation	C29H32N4O3S
Storage	3 years -20°C powder; 2 years -80°C in solvent;

## Images



## Product Description

Research Area: Others SMILES: CCS(=O)(=O)NC1=CC=C(C=C1)NC(=O)C2=C(/C3=CC=CC=C3)\NC4=CC=C(C=C4)CN5CCCCC5 Pathways: Cell Cycle/Checkpoint; Chromatin/Epigenetic Receptor: Aurora B (human); TbAUK1 Boling pt: Melting pt: Solubility: DMSO: 103 mg/mL (199.36 mM) Ethanol:

## References

1. Hauf S, et al. J Cell Biol, 2003, 161(2), 281-294.

Note: This product is for in vitro research use only