

Hexestrol


Catalog No: #S6533

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Hexestrol
Brief Description	Inhibitors
Purification	/
Target Name	Estrogen/progestogen Receptor agonist
Calculated MW	270.37
Formulation	C18H22O2
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area:Endocrine systemSMILES:CC[C@@H]([C@@H](CC)C1=CC=C(O)C=C1)C1=CC=C(O)C=C1

|r,c:11,19,t:6,8,14,16|Pathways:Endocrinology/HormonesReceptor:ERα; ERβBoiling pt:Melting pt:Solubility:DMSO:42 mg/mL (155.34 mM)Ethanol:42 mg/mL (155.34 mM)Water:

Note: This product is for in vitro research use only