

Histamine 2HCl


Catalog No: #S6534

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Histamine 2HCl
Brief Description	Inhibitors
Purification	99.00%
Target Name	Histamine Receptor agonist
Calculated MW	184.07
Formulation	C ₅ H ₉ N ₃ .2HCl
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Nervous system
 SMILES: Cl.C1.NCCC1=CNC=N1 [c:6,t:3]
 Pathways: GPCR/G Protein; Immunology/Inflammation;
 Neuroscience
 Receptor: H1 receptor; H2 receptor
 Boiling pt: Melting pt: Solubility: DMSO: 3 mg/mL (16.29 mM) Ethanol:

References

1. Hornyak SC, et al. Inflammation, 2003, 27(5), 317-327.

Note: This product is for in vitro research use only