

# Histamine Phosphate


Catalog No: #S6535

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Histamine Phosphate
Brief Description	Inhibitors
Purification	/
Target Name	Histamine Receptor agonist
Calculated MW	307.14
Formulation	C <sub>5</sub> H <sub>9</sub> N.2H <sub>3</sub> O <sub>4</sub> P
Storage	3 years -20°C powder; 2 years -80°C in solvent;

## Images



## Product Description

Research Area: Nervous system  
 SMILES: OP(O)(O)=O.OP(O)(O)=O.NCCC1=CNC=N1 | c:14,t:11 | Pathways: Neuroscience  
 Receptor: Histamine H1 receptor, Histamine H2 receptor  
 Boiling pt: Melting pt: Solubility: DMSO:

Note: This product is for in vitro research use only