

Ibandronate sodium monohydrate


Catalog No: #S6539

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Ibandronate sodium monohydrate
Brief Description	Inhibitors
Purification	/
Target Name	Others
Calculated MW	360.23
Formulation	C ₉ H ₂₃ NO ₇ P ₂ .Na.H ₂ O
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Other
 SMILES: O.[Na].CCCCCN(C)CCC(O)(P(O)(O)=O)P(O)(O)=O
 Pathways: Others
 Receptor: Others
 Boiling pt: Melting
 pt: Solubility: DMSO:

Note: This product is for in vitro research use only