

## Imidapril HCl

Catalog No: #S6542


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	Imidapril HCl
Brief Description	Inhibitors
Purification	/
Target Name	RAAS inhibitor
Calculated MW	441.91
Formulation	C20H27N3O6.HCl
Storage	3 years -20°C powder;2 years -80°C in solvent;

## Images



## Product Description

Research Area:Endocrine

systemSMILES:CCOC(=O)[C@H](CCC1=CC=CC=C1)N[C@@H](C)C(=O)N2[C@@H](CN(C2=O)C)C(=O)O.ClPathways:Endocrinology/

HormonesReceptor:ACEBoiling pt:Melting pt:Solubility:DMSO:88 mg/mL (199.13 mM)Ethanol:29 mg/mL (65.62 mM)Water:54 mg/mL (122.19

mM)Appearance:Remark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

Note: This product is for in vitro research use only