

## Isoconazole nitrate

Catalog No: #S6550


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	Isoconazole nitrate
Brief Description	Inhibitors
Purification	98.00%
Target Name	Others
Calculated MW	479.14
Formulation	C <sub>18</sub> H <sub>14</sub> Cl <sub>4</sub> N <sub>2</sub> O <sub>2</sub> HNO <sub>3</sub>
Storage	3 years -20°C powder; 2 years -80°C in solvent;

## Images



## Product Description

## Research

Area: Other  
 SMILES: O=[N+]([O-])=O.C1=CC(Cl)=C(C=C1)C(CN1C=CN=C1)OCC1=C(Cl)C=CC=C1Cl  
 Pathways: Others  
 Receptor: Others  
 Boiling pt: Melting pt: Solubility: Ethanol: 4 mg/mL (8.34 mM); DMSO: 41 mg/mL (85.56 mM)  
 Appearance: Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several months.

## References

1. Veraldi S. *Mycoses*. 2013 May; 56 Suppl 1: 41-3.

Note: This product is for in vitro research use only