

Isosorbide


Catalog No: #S6551

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Isosorbide
Brief Description	Inhibitors
Purification	/
Target Name	Others
Calculated MW	146.14
Formulation	C6H10O4
Storage	3 years -20°C powder;2 years -80°C in solvent;

Images



Product Description

Research Area:OtherSMILES:[H][C@]12OC[C@@H](O)[C@@]1([H])OC[C@@H]2O |r|Pathways:OthersReceptor:OthersBoiling pt:Melting pt:Solubility:DMSO:29 mg/mL (198.43 mM)Ethanol:29 mg/mL (198.43 mM)Water:29 mg/mL (198.43 mM)Appearance:Remark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

Note: This product is for in vitro research use only