

K-Ras(G12C) inhibitor 12

Catalog No: #S6555


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	K-Ras(G12C) inhibitor 12
Brief Description	Inhibitors
Purification	/
Target Name	Raf inhibitor
Calculated MW	449.67
Formulation	C15H17ClIN3O3
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Cancer SMILES: OC1=CC(Cl)=C(I)C=C1NCC(=O)N1CCN(CC1)C(=O)C=C |c:7,t:1,4|Pathways: Cell

Cycle/Checkpoint/Receptor: K-Ras(G12C) Boling pt: Melting pt: Solubility: DMSO: 27 mg/mL (60.04 mM) Ethanol:

Note: This product is for in vitro research use only