

Laquinimod

Catalog No: #S6561


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Laquinimod
Brief Description	Inhibitors
Purification	98.84%
Target Name	Others
Calculated MW	356.8
Formulation	C19H17ClN2O3
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Other SMILES: CCN(C(=O)C1=C(O)C2=C(C=CC=C2Cl)N(C)C1=O)C1=CC=CC=C1

|c:5,10,12,23,25,t:8,21|Pathways: Others Receptor: Others Boling pt: Melting pt: Solubility: DMSO: 61 mg/mL (170.96 mM) Ethanol: 1 mg/mL (2.8 mM) Water:

References

1. Zou LP, et al. Neuropharmacology, 2002, 42(5), 731-739.

Note: This product is for in vitro research use only