

Levobetaxolol HCl

Catalog No: #S6565


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Levobetaxolol HCl
Brief Description	Inhibitors
Purification	/
Target Name	Adrenergic Receptor antagonist
Calculated MW	343.89
Formulation	C18H29NO3.HCl
Storage	3 years -20°C powder;2 years -80°C in solvent;

Images



Product Description

Research Area:Nervous systemSMILES:CC(C)NC[C@@H](COC1=CC=C(C=C1)CCOCC2CC2)O.CIPathways:GPCR/G

ProteinReceptor:β1-adrenergic receptor , β2-adrenergic receptorBoiling pt:Melting pt:Solubility:DMSO:69 mg/mL (200.64 mM)Ethanol:69 mg/mL

(200.64 mM)Water:69 mg/mL (200.64 mM)Appearance:Remark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

Note: This product is for in vitro research use only