

Mucin 13 Polyclonal Antibody

Catalog No: #41189

Package Size: #41189-1 50ul #41189-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Mucin 13 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB IHC IF ELISA
Species Reactivity	Hu
Specificity	Mucin 13 Polyclonal Antibody detects endogenous levels of Mucin 13 protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the C-terminal region of human Mucin 13.
Target Name	Mucin 13
Other Names	MUC13; DRCC1; RECC; Mucin-13; MUC-13; Down-regulated in colon cancer 1
Accession No.	Swiss-Prot: Q9H3R2NCBI Gene ID: 56667
Uniprot	Q9H3R2
GeneID	56667;
SDS-PAGE MW	50kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details

Western Blot: 1/500 - 1/2000.


Immunohistochemistry: 1/100 - 1/300.

Immunofluorescence: 1/200 - 1/1000.


ELISA: 1/20000.

Not yet tested in other applications.

Images



Western Blot analysis of COLO205 cells using Mucin 13 Polyclonal Antibody



Western Blot analysis of L929 cells using Mucin 13 Polyclonal Antibody

Note: This product is for in vitro research use only