


## Description

Product Name	ML133 HCl
Brief Description	Inhibitors
Purification	/
Target Name	Potassium Channel inhibitor
Calculated MW	313.82
Formulation	C19H19NO.HCl
Storage	3 years -20°C powder;2 years -80°C in solvent;

## Images



## Product Description

Research Area:Cardiovascular SystemSMILES:Cl.COC1=CC=C(CNCC2=CC=CC3=C2C=CC=C3)C=C1

[c:11,13,16,18,21,t:2,4,9]Pathways:Membrane transporter/Ion channelReceptor:Kir2.1Boiling pt:Melting pt:Solubility:DMSO:63 mg/mL (200.75 mM)Ethanol:6 mg/mL (19.11 mM)Water:

Note: This product is for in vitro research use only