

Morin Hydrate


Catalog No: #S6596

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Morin Hydrate
Brief Description	Inhibitors
Purification	/
Target Name	Others
Calculated MW	320.25
Formulation	C ₁₅ H ₁₀ O ₇ ·H ₂ O
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Other SMILES: O.CC1=CC=C(C(O)=C1)C1=C(O)C(=O)C2=C(O)C=C(O)C=C2O1

|c:3,6,9,14,20,t:1,17|Pathways: Others Receptor: Others Boiling pt: Melting pt: Solubility: DMSO: 64 mg/mL (199.84 mM) Ethanol: 2 mg/mL warmed (6.24 mM) Water:

Note: This product is for in vitro research use only