

# N6-methyladenosine (m6A)


Catalog No: #S6599

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	N6-methyladenosine (m6A)
Brief Description	Inhibitors
Purification	99.00%
Target Name	Others
Calculated MW	281.27
Formulation	C11H15N5O4
Storage	3 years -20°C powder; 2 years -80°C in solvent;

## Images



## Product Description

Research Area: Other SMILES: CNC1=NC=NC2=C1N=CN2[C@@H]1O[C@H](CO)[C@@H](O)[C@H]1O

[r,c:4,6,9,t:2] Pathways: Others Receptor: Others Boiling pt: Melting pt: Solubility: DMSO: 56 mg/mL (199.09 mM) Ethanol: 2 mg/mL (7.11 mM) Water: 14 mg/mL (49.77 mM) Appearance: Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

## References

1. Aloni Y, et al. J Virol, 1979, 32(1), 52-60.

Note: This product is for in vitro research use only