

## Naftopidil DiHCl


Catalog No: #S6600

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Naftopidil DiHCl
Brief Description	Inhibitors
Purification	/
Target Name	5-HT Receptor antagonist; Adrenergic Receptor antagonist
Calculated MW	465.41
Formulation	C <sub>24</sub> H <sub>28</sub> N <sub>2</sub> O <sub>3</sub> ·2HCl
Storage	3 years -20°C powder; 2 years -80°C in solvent;

## Images



## Product Description

Research Area: Nervous system SMILES: Cl.Cl.COC1=CC=CC=C1N1CCN(CC(O)COC2=C3C=CC=CC3=CC=C2)CC1

[c:4,6,18,20,22,25,27,t:2] Pathways: GPCR/G Protein; Neuroscience Receptor: 5-HT1A; α1-adrenergic receptor Boling pt: Melting pt: Solubility: DMSO: 41 mg/mL (88.09 mM) Ethanol: 2 mg/mL (4.29 mM) Water:

Note: This product is for in vitro research use only