

## Nicaraven

Catalog No: #S6606


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	Nicaraven
Brief Description	Inhibitors
Purification	99.00%
Target Name	Others
Calculated MW	284.31
Formulation	C15H16N4O2
Storage	3 years -20°C powder;2 years -80°C in solvent;

## Images



## Product Description

Research Area:OtherSMILES:CC(CNC(=O)C1=CN=CC=C1)NC(=O)C1=CN=CC=C1 |c:8,10,18,20,t:6,16|Pathways:OthersReceptor:OthersBoiling pt:Melting pt:Solubility:DMSO:57 mg/mL (200.48 mM)Ethanol:57 mg/mL (200.48 mM)Water:57 mg/mL (200.48 mM)Appearance:Remark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

## References

1. Zingarelli B, et al. Shock, 2000, 13(2), 126-134.

Note: This product is for in vitro research use only