

## Nitenpyram

Catalog No: #S6607


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	Nitenpyram
Brief Description	Inhibitors
Purification	/
Target Name	AChR agonist
Calculated MW	270.72
Formulation	C11H15ClN4O2
Storage	3 years -20°C powder; 2 years -80°C in solvent;

## Images



## Product Description

Research Area: Nervous system  
 SMILES: CCN(CC1=CN=C(C=C1)Cl)/C(=C/[N+](=O)[O-])/N  
 Pathways: Neuroscience  
 Receptor: nAChR  
 Boiling pt: Melting pt: Solubility: DMSO: 54 mg/mL (199.46 mM) Ethanol: 54 mg/mL (199.46 mM) Water: 54 mg/mL warmed (199.46 mM)  
 Appearance: Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several months.

Note: This product is for in vitro research use only