

## Ondansetron


Catalog No: #S6615

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Ondansetron
Brief Description	Inhibitors
Purification	/
Target Name	5-HT Receptor antagonist
Calculated MW	293.36
Formulation	C18H19N3O
Storage	3 years -20°C powder; 2 years -80°C in solvent;

## Images



## Product Description

Research Area: Nervous system  
SMILES: CC1=NC=CN1CC2CCC3=C(C2=O)C4=CC=CC=C4N3C  
Pathways: Neuroscience  
Receptor: 5-HT<sub>3</sub> receptor  
Boiling pt: Melting pt: Solubility: DMSO: 0.1 mg/mL (0.34 mM) Ethanol:

Note: This product is for in vitro research use only