

## Ondansetron HCl

Catalog No: #S6616


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	Ondansetron HCl
Brief Description	Inhibitors
Purification	/
Target Name	5-HT Receptor antagonist
Calculated MW	329.82
Formulation	C18H19N3O.HCl
Storage	3 years -20°C powder;2 years -80°C in solvent;

## Images



## Product Description

Research Area:Nervous systemSMILES:Cl.CN1C2=CC=CC=C2C2=C1CCC(CN1C=CN=C1C)C2=O |c:4,6,9,17,19,t:2|Pathways:GPCR/G Protein;  
NeuroscienceReceptor:5-HT3Boiling pt:Melting pt:Solubility:DMSO:66 mg/mL (200.1 mM)Ethanol:10 mg/mL (30.31 mM)Water:24 mg/mL (72.76  
mM)Appearance:Remark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution  
can be stored below -20 °C for several month.

Note: This product is for in vitro research use only