

## Optovin

Catalog No: #S6617


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	Optovin
Brief Description	Inhibitors
Purification	99.00%
Target Name	TRP/TRPV Channel activator
Calculated MW	315.41
Formulation	C15H13N3OS2
Storage	3 years -20°C powder; 2 years -80°C in solvent;

## Images



## Product Description

Research Area: Other SMILES: CC1=CC(=C(N1C2=CN=CC=C2)C)/C=C\3/C(=O)NC(=S)S3

Pathways: Membrane transporter/Ion channel Receptor: TRPA1 Boling pt: Melting pt: Solubility: DMSO: 63 mg/mL warmed (199.74 mM) Ethanol:

## References

1. Kokel D, et al. Nat Chem Biol. 2013, 9(4), 257-263.

Note: This product is for in vitro research use only