

## Oxytetracycline Dihydrate

Catalog No: #S6624


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	Oxytetracycline Dihydrate
Brief Description	Inhibitors
Purification	/
Target Name	Others
Calculated MW	496.46
Formulation	C <sub>22</sub> H <sub>24</sub> N <sub>2</sub> O <sub>9</sub> ·2H <sub>2</sub> O
Storage	3 years -20°C powder; 2 years -80°C in solvent;

## Images



## Product Description

## Research

Area: OtherSMILES: O.O.[H][C@@]12[C@@H](O)[C@@]3([H])C(C(=O)C4=C(O)C=CC=C4[C@@]3(C)O)=C(O)[C@]1(O)C(=O)C(C(N)=O)=C(O)[C@H]2N(C)C |r,c:9,12,14,t:20,31|Pathways: OthersReceptor: OthersBoiling pt: Melting pt: Solubility: DMSO: 5 mg/mL (10.07 mM) Ethanol:

Note: This product is for in vitro research use only