

Ozagrel hydrochloride

Catalog No: #S6625


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Ozagrel hydrochloride
Brief Description	Inhibitors
Purification	98.00%
Target Name	Thrombin inhibitor
Calculated MW	264.71
Formulation	C ₁₃ H ₁₂ N ₂ O ₂ HCl
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Metabolism system
 SMILES: C1=CC(=CC=C1CN2C=CN=C2)/C=C/C(=O)O
 Pathways: Others
 Receptor: Thromboxane A₂ synthetase
 Boiling pt: Melting pt: Solubility: Ethanol: 6 mg/mL (22.66 mM); Water: 53 mg/mL (200.21 mM); DMSO: 53 mg/mL (200.21 mM)
 Appearance: Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several months.

References

1. Ishitsuka Y, et al. J Pharmacol Sci, 2009, 111(2), 211-215.

Note: This product is for in vitro research use only