

Phosphoramidon Disodium Salt

Catalog No: #S6627


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Phosphoramidon Disodium Salt
Brief Description	Inhibitors
Purification	/
Target Name	Others
Calculated MW	588.48
Formulation	C ₂₃ H ₃₄ N ₃ Na ₂ O ₁₀ P
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research

Area: Other SMILES: [Na][Na].CC(C)C[C@H](NP(O)(=O)O)[C@@H]1O[C@H](C)[C@H](O)[C@@H](O)[C@H]1O)C(=O)N[C@@H](CC1=CNC2=CC=CC=C12)C(O)=O |r,c:32,t:27,30,34|Pathways: Others Receptor: metalloendopeptidase Boling pt: Melting pt: Solubility: DMSO: 100 mg/mL (169.92 mM) Ethanol:

Note: This product is for in vitro research use only