

Quercetin Dihydrate

Catalog No: #S6630


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Quercetin Dihydrate
Brief Description	Inhibitors
Purification	99.00%
Target Name	Others
Calculated MW	338.27
Formulation	C ₁₅ H ₁₀ O ₇ ·2H ₂ O
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Cardiovascular

system SMILES: O.O.OC1=CC(O)=C2C(OC(=C(O)C2=O)C2=CC(O)=C(O)C=C2)=C1 Pathways: Others Receptor: Others Boling pt: Melting

pt: Solubility: DMSO: 68 mg/mL (201.02 mM) Ethanol: 21 mg/mL (62.08 mM) Water:

References

1. Dajas F, J Ethnopharmacol, 2012, 143(2), 383-396.

Note: This product is for in vitro research use only