

Resminostat


Catalog No: #S6637

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Resminostat
Brief Description	Inhibitors
Purification	98.00%
Target Name	HDAC inhibitor
Calculated MW	349.4
Formulation	C16H19N3O4S
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Cancer
 SMILES: CN(C)CC1=CC=C(C=C1)S(=O)(=O)N2C=CC(=C2)/C=C/C(=O)NO
 Pathways: Chromatin/Epigenetic; DNA Damage/DNA Repair; NF-kb Receptor: HDAC1; HDAC3; HDAC6
 Boling pt: Melting pt: Solubility: DMSO: 70 mg/mL (200.34 mM) Ethanol: 70 mg/mL (200.34 mM) Water:

Note: This product is for in vitro research use only