

Retapamulin

Catalog No: #S6638


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Retapamulin
Brief Description	Inhibitors
Purification	98.00%
Target Name	ribosome inhibitor
Calculated MW	517.76
Formulation	C30H47NO4S
Storage	3 years -20°C powder;2 years -80°C in solvent;

Images



Product Description

Research

Area: OtherSMILES: [H][C@]12CC[C@]([H])(C[C@H](C1)SCC(=O)O[C@@H]1C[C@@](C)(C=C)[C@H](O)[C@H](C)[C@]34CCC(=O)[C@@]3([H])[C@@]1(C)[C@H](C)CC4)N2C |r|Pathways: Microbiology&virologyReceptor: ribosomeBoiling pt: Melting pt: Solubility: DMSO: 104 mg/mL (200.86 mM)Ethanol: 104 mg/mL (200.86 mM)Water:

References

1. Yan K, et al. Antimicrob Agents Chemother, 2006, 50(11), 3875-3881.

Note: This product is for in vitro research use only