

Rimantadine


Catalog No: #S6642

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Rimantadine
Brief Description	Inhibitors
Purification	/
Target Name	Influenza Virus inhibitor
Calculated MW	179.3
Formulation	C12H21N
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Immunology
SMILES: CC(N)C12CC3CC(CC(C3)C1)C2
Pathways: Microbiology & virology
Receptor: Influenza virus
Boiling pt: Melting pt: Solubility: DMSO: 33 mg/mL (184.04 mM) Ethanol: 33 mg/mL (184.04 mM) Water: 33 mg/mL (184.04 mM)
Appearance: Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several months.

Note: This product is for in vitro research use only