

Ro 31-8220 Mesylate

Catalog No: #S6643


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Ro 31-8220 Mesylate
Brief Description	Inhibitors
Purification	99.00%
Target Name	PKC inhibitor
Calculated MW	553.65
Formulation	C25H23N5O2SCH4O3S
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research

Area: Immunology SMILES: CS(O)(=O)=O.CN1C=C(C2=CC=CC=C2)C1=C(C(=O)NC1=O)C1=CN(CCCSC(N)=N)C2=CC=CC=C12 Pathways: Cytoskeletal Signaling Receptor: PKC α ; PKC β 1; PKC β 2; PKC γ ; PKC ϵ Boiling pt: Melting pt: Solubility: DMSO: 100 mg/mL (180.61 mM) Appearance: Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several months.

References

1. Wilkinson SE, et al. Biochem J. 1993, 294 (Pt 2), 335-337.

Note: This product is for in vitro research use only