

Rotigotine


Catalog No: #S6647

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Rotigotine
Brief Description	Inhibitors
Purification	99.59%
Target Name	Dopamine Receptor agonist
Calculated MW	315.47
Formulation	C19H25NOS
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Nervous system
 SMILES: CCCN(CCC1=CC=CS1)[C@H]1CCC2=C(C1)C=CC=C2O | r,c:8,15,19,21,t:6 | Pathways: GPCR/G Protein;
 Neuroscience Receptor: D2; D3
 Boling pt: Melting pt: Solubility: DMSO: 63 mg/mL (199.7 mM) Ethanol: 63 mg/mL (199.7 mM) Water:

References

1. Giladi N, et al. J Parkinsons Dis. 2016 Apr 2

Note: This product is for in vitro research use only