

Safinamide Mesylate

Catalog No: #S6651


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Safinamide Mesylate
Brief Description	Inhibitors
Purification	99.85%
Target Name	MAO inhibitor
Calculated MW	398.45
Formulation	C17H19FN2O2.CH4O3S
Storage	3 years -20°C powder;2 years -80°C in solvent;

Images



Product Description

Research Area:Nervous systemSMILES:CS(O)(=O)=O.C[C@H](NCC1=CC=C(OCC2=CC(F)=CC=C2)C=C1)C(N)=O

[r,c:17,19,22,t:8,10,14]Pathways:MetabolismReceptor:MAO-B , MAO-ABoling pt:Melting pt:Solubility:DMSO:80 mg/mL (200.77 mM);Ethsbol:13 mg/mL (32.62 mM);H2O:80 mg/mL (200.77 mM)Appearance:Remark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

Note: This product is for in vitro research use only