

## Sapogenins Glycosides

Catalog No: #S6654


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	Sapogenins Glycosides
Brief Description	Inhibitors
Purification	/
Target Name	Others
Calculated MW	388.95
Formulation	C18H28N2O3S.HCl
Storage	3 years -20°C powder;2 years -80°C in solvent;

## Images



## Product Description

## Research

Area: OtherSMILES: C/C=C(/C)\C(=O)O[C@H]1[C@@H]([C@@]2[C@@H](C[C@@]3C(=CCC4[C@]3(CCC5[C@@]4(CC[C@@H]([C@]5(C)CO)O[C@H]6[C@@H]([C@H]([C@@H]([C@H](O6)C(=O)O)O[C@H]7[C@@H]([C@H]([C@@H]([C@H](O7)CO)O)O)O)O[C@H]8[C@@H]([C@H]([C@@H]([C@H](O8)CO)O)O)O)C)C)C)C)O)CO)OC(=O)CPathways: OthersReceptor: OthersBoiling pt: Melting pt: Solubility: N/AAppearance: Remark: For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

Note: This product is for in vitro research use only