

Scopine


Catalog No: #S6662

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Scopine
Brief Description	Inhibitors
Purification	97.00%
Target Name	Adrenergic Receptor agonist
Calculated MW	155.19
Formulation	C8H13NO2
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Others SMILES: CN1[C@@H]2CC(C[C@H]1[C@H]3[C@@H]2O3)O Pathways: GPCR/G Protein Receptor: α 1-adrenergic receptor Boiling pt: Melting pt: Solubility: DMSO: 31 mg/mL (199.75 mM) Ethanol:

References

1. Wang X, et al. *Bioconjug Chem.* 2014 Nov 19;25(11):2046-54.

Note: This product is for in vitro research use only