

Scopine HCl


Catalog No: #S6663

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Scopine HCl
Brief Description	Inhibitors
Purification	/
Target Name	Adrenergic Receptor agonist
Calculated MW	191.66
Formulation	C8H13NO2.HCl
Storage	3 years -20°C powder;2 years -80°C in solvent;

Images



Product Description

Research Area:Nervous systemSMILES:CN1[C@@H]2CC(C[C@H]1[C@H]3[C@@H]2O3)O.ClPathways:NeuroscienceReceptor:α1-adrenergic receptor Boling pt:Melting pt:Solubility:DMSO:31 mg/mL (199.75 mM)Ethanol:

Note: This product is for in vitro research use only