

SecinH3

Catalog No: #S6664


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	SecinH3
Brief Description	Inhibitors
Purification	99.00%
Target Name	Others inhibitor
Calculated MW	460.5
Formulation	C ₂₄ H ₂₀ N ₄ O ₄ S
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Other SMILES: COC1=NN(C(=N1)C1=CC=C2OCOC2=C1)C1=CC=C(NC(=O)CSC2=CC=CC=C2)C=C1

|c:5,16,30,32,35,t:2,8,10,19,21,28|Pathways: Others Receptor: drosophila steppke; hCyh1; hCyh2; hCyh3; mCyh3 Boling pt: Melting

pt: Solubility: DMSO: 92 mg/mL (199.78 mM) Ethanol: 1 mg/mL warmed (2.17 mM) Water:

References

1. Hafner M, et al. Nature. 2006, 444(7121), 941-944.

Note: This product is for in vitro research use only