

## Silymarin

Catalog No: #S6670


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	Silymarin
Brief Description	Inhibitors
Purification	98.00%
Target Name	Others
Calculated MW	482.44
Formulation	C25H22O10
Storage	3 years -20°C powder;2 years -80°C in solvent;

## Images



## Product Description

## Research

Area:OtherSMILES:COC1=C(O)C=CC(=C1)C1OC2=C(OC1CO)C=C(C=C2)C1OC2=CC(O)=CC(O)=C2C(=O)C1OPathways:OthersReceptor:OthersB  
 oling pt:Melting pt:Solubility:Ethanol: 1 mg/mL (2.07 mM); DMSO: 97 mg/mL (201.06 mM)Appearance:YellowPowderRemark:For obtaining a higher  
 solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

## References

1. Feher J, et al. Curr Pharm Biotechnol, 2012, 13(1), 210-217.

Note: This product is for in vitro research use only