

Sodium Picosulfate

Catalog No: #S6675


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Sodium Picosulfate
Brief Description	Inhibitors
Purification	98.00%
Target Name	Others inhibitor
Calculated MW	481.41
Formulation	C ₁₈ H ₁₅ NO ₈ S ₂ .2Na
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area:OtherSMILES:[Na+].[Na+].[O-]S(=O)(=O)OC1=CC=C(C=C1)C(C1=CC=C(OS(=O)(=O)C=C1)C1=CC=CC=N1

|c:7,9,22,27,29,t:5,13,15,25|Pathways:OthersReceptor:OthersBoiling pt:Melting pt:Solubility:DMSO:96 mg/mL (199.41 mM)Ethanol:3 mg/mL (6.23 mM)Water:96 mg/mL (199.41 mM)Appearance:Remark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

References

1. Nishikawa J, et al. Arzneimittelforschung, 1981, 31(2), 321-325.

Note: This product is for in vitro research use only