

## Succinylcholine Chloride Dihydrate

Catalog No: #S6684


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	Succinylcholine Chloride Dihydrate
Brief Description	Inhibitors
Purification	/
Target Name	AChR agonist
Calculated MW	397.34
Formulation	C <sub>14</sub> H <sub>30</sub> N <sub>2</sub> O <sub>4</sub> ·2Cl·2H <sub>2</sub> O
Storage	3 years -20°C powder; 2 years -80°C in solvent;

## Images



## Product Description

Research Area: Nervous system  
 SMILES: O.O.[Cl-].[Cl-].C[N+](C)(C)CCOC(=O)CCC(=O)OCC[N+](C)(C)C  
 Pathways: Neuroscience  
 Receptor: nicotinic AChR  
 Boiling pt: Melting pt: Solubility: DMSO: 1 mg/mL (2.51 mM) Ethanol:

Note: This product is for in vitro research use only