

Suprofen

Catalog No: #S6687


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Suprofen
Brief Description	Inhibitors
Purification	98.00%
Target Name	COX inhibitor
Calculated MW	260.31
Formulation	C14H12O3S
Storage	3 years -20°C powder;2 years -80°C in solvent;

Images



Product Description

Research Area:OtherSMILES:CC(C(O)=O)C1=CC=C(C=C1)C(=O)C1=CC=CS1 |c:7,9,16,t:5,14|Pathways:Immunology/Inflammation;
NeuroscienceReceptor:COX-1; COX-2Boiling pt:Melting pt:Solubility:10 mM in DMSOAppearance:Remark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

References

1.Yeadon, A., et al., Suprofen. An overview of long-term safety. Pharmacology, 1983. 27 Suppl 1: p. 87-94.

Note: This product is for in vitro research use only