


## Description

Product Name	T0901317
Brief Description	Inhibitors
Purification	99.00%
Target Name	FXR agonist; Liver X Receptor agonist
Calculated MW	481.33
Formulation	C17H12F9NO3S
Storage	3 years -20°C powder; 2 years -80°C in solvent;

## Images



## Product Description

Research Area: Metabolism system SMILES: OC(C1=CC=C(C=C1)N(CC(F)(F)F)S(=O)(=O)C1=CC=CC=C1)(C(F)(F)F)C(F)(F)F

|c:4,6,20,22,t:2,18|Pathways: Metabolism; Others Receptor: FXR; LXR Boling pt: Melting pt: Solubility: DMSO: 96 mg/mL (199.44 mM) Ethanol: 96 mg/mL (199.44 mM) Water:

## References

1. Repa JJ, et al. Science, 2000, 289(5484), 1524-1529.

Note: This product is for in vitro research use only