

TAME

Catalog No: #S6693


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	TAME
Brief Description	Inhibitors
Purification	/
Target Name	APC inhibitor
Calculated MW	342.41
Formulation	C14H22N4O4S
Storage	3 years -20°C powder;2 years -80°C in solvent;

Images



Product Description

Research Area:CancerSMILES:COC(=O)[C@H](CCCNC(N)=N)NS(=O)(=O)C1=CC=C(C)C=C1 |r,c:21,t:16,18|Pathways:Cell

Cycle/CheckpointReceptor:APCBoling pt:Melting pt:Solubility:DMSO:69 mg/mL (201.51 mM)Ethanol:3 mg/mL (8.76 mM)Water:69 mg/mL (201.51 mM)Appearance:Remark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

Note: This product is for in vitro research use only