

Tamibarotene

Catalog No: #S6694


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Tamibarotene
Brief Description	Inhibitors
Purification	99.00%
Target Name	Retinoid Receptor agonist
Calculated MW	351.44
Formulation	C22H25NO3
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Metabolism

system SMILES: CC1(C)CCC(C)(C)C2=C1C=CC(NC(=O)C1=CC=C(C=C1)C(=O)O)=C2 Pathways: Metabolism Receptor: RAR α ; RAR β Boling pt: Melting pt: Solubility: Ethanol: 70 mg/mL (199.18 mM); DMSO: 70 mg/mL (199.18 mM) Appearance: Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several months.

References

1. Sanda T, et al. Leukemia, 2005, 19(6), 901-909.

Note: This product is for in vitro research use only