

## Tazarotene

Catalog No: #S6696


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	Tazarotene
Brief Description	Inhibitors
Purification	97.01%
Target Name	Retinoid Receptor inhibitor
Calculated MW	351.46
Formulation	C <sub>21</sub> H <sub>21</sub> NO <sub>2</sub> S
Storage	3 years -20°C powder; 2 years -80°C in solvent;

## Images



## Product Description

Research Area: Other  
 SMILES: CCOC(=O)C1=CN=C(C=C1)C#CC1=CC2=C(SCCC2(C)C)C=C1  
 Pathways: Metabolism  
 Receptor: RAR  
 Boiling pt: Melting pt: Solubility: Ethanol: 10 mg/mL (28.45 mM); DMSO: 10 mg/mL (28.45 mM)  
 Appearance: Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

## References

1. Yen A, et al. Mol Pharmacol, 2004, 66(6), 1727-1737.

Note: This product is for in vitro research use only