

Tenatoprazole

Catalog No: #S6700


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Tenatoprazole
Brief Description	Inhibitors
Purification	99.50%
Target Name	Proton pump inhibitor
Calculated MW	346.4
Formulation	C16H18N4O3S
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Cardiovascular System SMILES: COC1=CC=C2NC(=NC2=N1)S(=O)CC1=C(C)C(OC)=C(C)C=N1

[c:7,10,16,24,t:2,4,21] Pathways: Membrane transporter/Ion channel Receptor: Proton pump Boiling pt: Melting pt: Solubility: DMSO: 46 mg/mL (132.79 mM) Ethanol: 14 mg/mL warmed (40.41 mM) Water:

References

1. Shin JM, et al. Biochem Pharmacol, 2006, 71(6), 837-849.

Note: This product is for in vitro research use only